

Biohazardous Waste - Do's and Don'ts

Do	Don't
 <p data-bbox="521 453 808 552">Dispose of blunt tip pipettes in a biomedical waste container.</p>	 <p data-bbox="1208 453 1471 552">Don't pre-treat disposable pipettes in a graduated cylinder.</p>
 <p data-bbox="558 800 808 1037">Use properly labeled biomedical waste containers with the biohazard symbol or the words "Biohazardous Waste".</p>	 <p data-bbox="1224 800 1455 856">Don't use non-labeled containers.</p>
 <p data-bbox="532 1178 800 1241">Tie off biomedical bag when $\frac{3}{4}$ full.</p>	 <p data-bbox="1256 1178 1417 1241">Don't overfill container.</p>
 <p data-bbox="532 1530 800 1593">Keep lid closed except when adding waste.</p>	 <p data-bbox="1224 1530 1466 1562">Don't keep lid open.</p>

Do not dispose of biohazardous waste in the regular trash.